

Please read the information contained herein carefully.



FreshKon® CLEAR ADVANCE Multi-Purpose Solution

FreshKon® CLEAR ADVANCE Multi-Purpose Solution (MPS) cleans, rinses, removes protein, disinfects, lubricates and stores all soft (hydrophilic) contact lenses including cosmetic and Silicone Hydrogel contact lenses.

FreshKon® CLEAR ADVANCE Multi-Purpose Solution comes with an additional exclusive wetting agent - Sodium Hyaluronate, which helps to further enhance moisture and lens wear comfort with its humidifying effect.

Warnings

It is important to follow the instruction of your eye-care practitioner and the directions as indicated on the package and label for proper handling and lens care. Failure to do so may result in severe eye infection or loss of vision. If you experience any eye discomfort, irritation, pain and redness or blurred vision, remove the lenses and consult your eye-care practitioner IMMEDIATELY.

DO NOT MIX WITH OTHER FLUIDS EXCEPT AS DIRECTED.

Do not use concurrently with any eye medication except on medical advice.

Adverse Reaction

The following problems may occur: Stinging or burning sensation in the eyes, irritation, discomfort, foreign body sensation, excessive tearing, redness, blurred vision, halos or ghosting around images, sensitivity to light or dry eyes.

If you experience any of the above, remove the lenses and consult your eye-care practitioner IMMEDIATELY.

Contraindication

If you are allergic to any ingredients in FreshKon® CLEAR ADVANCE Multi-Purpose Solution, do not use this solution.

Precautions

- Do not purchase or use FreshKon® CLEAR ADVANCE Multi-Purpose Solution if bottle is opened or tempered/damaged.
- Do not reuse FreshKon® CLEAR ADVANCE Multi-Purpose Solution.
- Do not use with any heat (thermal) disinfecting method.
- To avoid contamination, do not touch the dropper tip of the container to any surface.
- Keep the container tightly capped when not in use.
- Store at room temperature (15°C - 30°C) and away from direct sunlight.
- Discard any remaining solution 90 days after first opening the bottle.
- Use before the expiry date.
- Keep out of reach of children.
- Do not freeze.

Instructions

- Wash hands with oil-free soap and dry with lint-free towel before handling contact lenses.
- Apply few drops of FreshKon® CLEAR ADVANCE Multi-Purpose Solution on lens and gently rub for 20 seconds.
- Rinse lens thoroughly with FreshKon® CLEAR ADVANCE Multi-Purpose Solution and soak for at least 6 hours before use.
- Rinse soaked lens with FreshKon® CLEAR ADVANCE Multi-Purpose Solution before wearing.
- Always empty lens case, clean with fresh solution and air dry before using.

Recommendation

Change your lens case and accessories at least once every 3 months.