

Instruction For Use



FreshKon® CLEAR Multi-Purpose Solution is a contact lens disinfection solution that cleans, rinses, removes protein, disinfects, lubricates and stores all soft (hydrophilic) contact lenses including cosmetic and Silicone Hydrogel contact lenses.

FreshKon® CLEAR is a Multi-Purpose Solution that has ionic molecules which combine with protein's charged molecules and enhance the release of protein deposits from the lens. HPMC forms a lubricating shield on the lens, keeps in moisture and lubricates the lens to give long lasting comfort to your eyes.

Contents

A sterile, aqueous, isotonic solution containing Poloxamer, Edetate Disodium (EDTA), Hypromellose (HPMC), Glycerine, Boric Acid, Borax, Sodium Chloride, Potassium Chloride and Polyhexanide (PHMB).

Warnings

It is important to follow the instructions of your eye-care practitioner and the directions as indicated on the package and label for proper handling and lens care. Failure to do so may result in severe eye infection or loss of vision. If you experience any eye discomfort, irritation, pain and redness or blurred vision, remove the lenses and consult your eye-care practitioner IMMEDIATELY.

DO NOT MIX WITH OTHER FLUIDS EXCEPT AS DIRECTED.

Do not use concurrently with any eye medication except on medical advice.

Adverse Reaction

The following problems may occur: Stinging or burning sensation in the eyes, irritation, discomfort, foreign body sensation, excessive tearing, redness, blurred vision, halos or ghosting around images, sensitivity to light or dry eyes.

If you experience any of the above, remove the lenses and consult your eye-care practitioner IMMEDIATELY.

Contraindication

If you are allergic to any ingredients in FreshKon® CLEAR Multi-Purpose Solution, do not use this solution.

Precautions

- Do not purchase or use FreshKon® CLEAR Multi-Purpose Solution if bottle is opened or tempered/damaged.
- Do not reuse or "top-off" FreshKon® CLEAR Multi-Purpose Solution.
- Do not use with any heat (thermal) disinfecting method.
- To avoid contamination, do not touch the dropper tip of the bottle to any surface.
- Keep the bottle tightly capped when not in use.
- Store at room temperature (15°C - 30°C) and away from direct sunlight.
- Discard any remaining solution 90 days after first opening the bottle.
- Use before the expiry date.
- Do not expose or rinse lenses or lens case with any water or non-sterile solution.

- Keep out of reach of children.
- Do not freeze.

Instructions

- Wash hands with oil-free soap and dry with lint-free towel before handling contact lenses.
- Apply few drops of FreshKon® CLEAR Multi-Purpose Solution on lens and gently rub for 20 seconds.
- Rinse lens thoroughly with FreshKon® CLEAR Multi-Purpose Solution and soak for at least 6 hours before next use.
- Rinse soaked lens with FreshKon® CLEAR Multi-Purpose Solution before wearing.
- Always empty lens case, clean with fresh solution and air dry before using.
- Change solution in lens case everyday even contact lens is not in use.
- It is advisable to consult an eye care practitioner before changing your lens care system.

Recommendation

Change your lens case and accessories at least once every 3 months.



Brand Owner
Oculus Private Limited
10 Kaki Bukit Avenue 1 #06-08
Singapore 417942
www.freshkon.com

Last revision July 2020

SFCM-IFUK(0)